

# NBIA 5.1 release announcement

Dear NBIA Users,

We're pleased to announce the release of version 5.1 of the National Biomedical Imaging Archive (NBIA), available at <http://imaging.nci.nih.gov>. This release significantly improves end user functionality and back end administration as explained below.

## Ultrasound Search

The NBIA Portal now has ultrasound specific searching capabilities.

Number of Frames (US)	<input checked="" type="checkbox"/> Single Frame <input checked="" type="checkbox"/> Multi-frame (Cine Loop)	
Ultrasound Color Data Present (US)	<input checked="" type="checkbox"/> Yes <input checked="" type="checkbox"/> No	
Ultrasound Image Mode (US)	Available	Selected
	<div><input type="checkbox"/> 2D Imaging</div> <div><input type="checkbox"/> Color Doppler</div> <div><input type="checkbox"/> M-Mode</div> <div><input type="checkbox"/> PW Doppler</div>	<div><div>SELECT ALL</div><div>REMOVE ALL</div></div>

[nci.nih.gov/pages/editpage.action?pageId=96109785#](https://wiki.nci.nih.gov/pages/editpage.action?pageId=96109785#)

[https://wiki.](https://wiki.nci.nih.gov/pages/editpage.action?pageId=96109785#)

The Lookup Search Terms in NBIA Simple Search retrieves the search terms for Ultrasound related criteria from remote nodes.

Mozilla Firefox

https://imaging-stage.nci.nih.gov/ncia/availableSearchTerms.jsf

NCICBIT-DEV

Modalities	CR, CT, DX, HC, HISTOPATHOLOGY, MR, PR, PT, RTDOSE, RTPLAN, RTSTRUCT, SC, SR, US, XA
Anatomic Sites	ABDOMEN, BRAIN, BREAST, CHEST, HEAD, LUNG, NECK, NOT SPECIFIED, RECTUM, THORAX_1_ROUTINE
Collections	CT Colonography, LIDC, NCRI, Phantom, Phantom-FDA, RIDER Breast MRI, RIDER Lung CT, RIDER Lung PET CT, RIDER Neuro MRI, RIDER Phantom MRI, RIDER Pilot, RoswellStrong, TCGA, VASARI, Virtual Colonoscopy
Convolution Kernels	7.200000 mm, B, B19f, B20f, B20s, B30f, B30s, B31f, B31s, B40s, B41f, B41s, B45f, B50f, B70f, B70s, B80f, B80s, BONE, C, D, DETAIL, EXPERIMENTAL7, FC01, FC02, FC03, FC10, FC11, FC13, FC50, FL01, FL02, H30s, LUNG, NULL, PET AC, Rad., SOFT, STANDARD, T20s, T80s, UA, UB, X-Y-Z Gaussian F, XYZ G5.00, XYZ G7.00, XYZ Gauss4.00, XYZ Gauss5.00, XYZ Gauss7.00, hanning, rectangle
Equipment	CMS, Inc., CPS, DeJarnette Research Systems, FUJI PHOTO FILM Co., Ltd., GE, GE MEDICAL SYSTEMS, GEMS, General Electric, KODAK, Multiple, N/A, Normalized to NAVM Cr, Philips, Philips Medical Systems, Philips Medical Systems (Netherlands), Philips Medical Systems Inc, PowerDicom, Radiology Research, SIEMENS, TOSHIBA, VARIAN Medical Systems, "GE Healthcare"
Ultrasound Image Mode	2D Imaging, M-Mode, PW Doppler, Color Doppler

NCICBIT-QA

Modalities	CR, CT, DX, HC, HISTOPATHOLOGY, MR, NM, OT, PR, PT, RTDOSE, RTSTRUCT, SC, SR, US
Anatomic Sites	ABDOMEN, BRAIN, BREAST, CHEST, HAND, HEAD, KNEE, LUNG, NECK, NOT SPECIFIED, PELVIS, RECTUM, THORAX_1_ROUTINE
Collections	Anderson, CT Colonography, Head-Neck Cetuximab, LIDC, NCRI, Phantom, Phantom-FDA, Project, Project_Sep, RIDER Breast MRI, RIDER Lung CT, RIDER Lung PET CT, RIDER Neuro MRI, RIDER Phantom MRI, RIDER Pilot, RoswellStrong, TCGA, VASARI, Virtual Colonoscopy
Convolution Kernels	7.200000 mm, B, B20f, B20s, B30f, B30s, B31f, B31s, B40s, B41f, B41s, B45f, B50f, B70f, B70s, B80f, B80s, BONE, C, D, DETAIL, EXPERIMENTAL7, FC01, FC02, FC03, FC10, FC11, FC13, FC50, FL01, FL02, H30s, LUNG, NULL, PET AC, Rad., SOFT, STANDARD, T20s, T80s, UA, UB, X-Y-Z Gaussian F, XYZ G5.00, XYZ G7.00, XYZ Gauss4.00, XYZ Gauss5.00

## QC Tool Improvement

The QC Tool now has a check all button to quickly select all items.

**QC Tool**

Search Controls

**SEARCH CRITERIA**

QC Status: ☒ Not Yet Reviewed ☒ Visible ☐ Not Visible ☐ To Be Deleted

Collection(s):

Patient ID:

Most recent submission date: From  To

**SUBMIT**

QC Images by Series

☒ CHECK ALL ☐ UNCHECK ALL

Results Per Page: 10

846 series found, displaying 10 records, from 31 to 40. Page 4 / 85

<input type="checkbox"/>	Submission date	Collection/Site	Patient	Study	Series	Visibility	Viewer
<input checked="" type="checkbox"/>	08/04/2006 00:00:00 AM	IDRI/IDRI	1.3.6.1.4.1.9328.50.3.0125	1.3.6.1.4.1.9328.50.3.38316	1.3.6.1.4.1.9328.50.3.38317	Visible	<a href="#">View</a>
<input checked="" type="checkbox"/>	08/04/2006 00:00:00 AM	IDRI/IDRI	1.3.6.1.4.1.9328.50.3.0127	1.3.6.1.4.1.9328.50.3.39182	1.3.6.1.4.1.9328.50.3.39183	Visible	<a href="#">View</a>
<input checked="" type="checkbox"/>	08/04/2006 00:00:00 AM	IDRI/IDRI	1.3.6.1.4.1.9328.50.3.0128	1.3.6.1.4.1.9328.50.3.39696	1.3.6.1.4.1.9328.50.3.39697	Visible	<a href="#">View</a>
<input checked="" type="checkbox"/>	08/04/2006 00:00:00 AM	IDRI/IDRI	1.3.6.1.4.1.9328.50.3.0129	1.3.6.1.4.1.9328.50.3.40151	1.3.6.1.4.1.9328.50.3.40152	Visible	<a href="#">View</a>
<input checked="" type="checkbox"/>	08/04/2006 00:00:00 AM	IDRI/IDRI	1.3.6.1.4.1.9328.50.3.0130	1.3.6.1.4.1.9328.50.3.40579	1.3.6.1.4.1.9328.50.3.40580	Visible	<a href="#">View</a>
<input checked="" type="checkbox"/>	08/04/2006 00:00:00 AM	IDRI/IDRI	1.3.6.1.4.1.9328.50.3.0131	1.3.6.1.4.1.9328.50.3.41039	1.3.6.1.4.1.9328.50.3.41040	Visible	<a href="#">View</a>
<input checked="" type="checkbox"/>	08/04/2006 00:00:00 AM	IDRI/IDRI	1.3.6.1.4.1.9328.50.3.0132	1.3.6.1.4.1.9328.50.3.41510	1.3.6.1.4.1.9328.50.3.41511	Visible	<a href="#">View</a>
<input checked="" type="checkbox"/>	08/04/2006 00:00:00 AM	IDRI/IDRI	1.3.6.1.4.1.9328.50.3.0133	1.3.6.1.4.1.9328.50.3.42013	1.3.6.1.4.1.9328.50.3.42014	Visible	<a href="#">View</a>

Improved Download Manager

Through collaborative work with the members of the Imaging Workspace, the reliability of Download Manager has been improved. If network errors occur during large downloads, Download Manager will retry the files a configurable number of times.

Improved Large Collection Deletion

The deletion of large collections is now performed asynchronously with a message driven bean. This allows the user to perform the deletion and go on to other tasks as the deletion is performed in the background. When the deletion is complete, the user receives an email detailing the aspects of the deletion.

Modality Information in Web Search

The text of the modality is displayed when it is in focus.

Image Modality(ies)

Available

☐ CR  
☐ CT  
☐ DX  
☐ HC  
☒ Digital Radiography  
☐ MR  
☐ NM  
☐ CT

SELECT ALL

REMOVE ALL

Selected\*

\*Return cases that include

☐ ALL
☒ ANY of the selected modalities.

Improved Security

NBIA has undergone extensive security scanning and is now protected for various forms of attack such as Cross Site Scripting and Cross Site Request Forgery,

Clinical Trials Processor (CTP) Update

NBIA has been integrated with the latest version of CTP

AntHill Pro 3

The NBIA can now be deployed with AntHill Pro 3, the tool used by CBIIT systems support.

Section 508 Compliance

The NBIA has improved compliance with Section 508 for people with disabilities.

## Technology Infrastructure Updates

The NBIA portal has been upgraded to run on JBoss 5.1, MySQL 5.1, and caGrid 1.4.

### For more information:

Detailed release notes can be found on [SVN](#)

For groups interested in installing NBIA 5.1 locally, the installation guide is located here [NBIA 5.1 Installation Guide](#)

As always, support is available from [Application Support](#).