

G-DOC Plus Help

How to perform gene expression based survival analysis

Innovation Center for Biomedical Informatics Georgetown University

Overview

- Login
- Navigation
- How to create patient groups
- Different ways of performing survival analysis on gene expression data (name of tool: Gene expression KM plot)
 - Going from Group Comparison to Gene expression KM
 - Create groups, perform Gene expression KM
 - Going to Gene expression KM directly run on all samples

Login



The Innovation Center for Biomedical Informatics (ICBI) Lombardi Comprehensive Cancer Center Thu Jan 22, 2015

kb472

•••••

Log In

register now | forgot password



Welcome to GDOC Plus Beta!

The Georgetown Database of Cancer Plus other diseases (G-DOC Plus) is a precision medicine platform containing molecular and clinical data from thousands of patients and cell lines, along with tools for analysis and data visualization. The platform enables the integrative analysis of multiple data types to understand disease

Precision Medicine

Translational research

Population genetics

Understanding Data in G-DOC Plus

It all begins with a study...

All data in G-DOC Plus derives from studies on topics such as breast cancer, wound healing, or even 1,000 Genomes. Each study may contain clinical and/or biospecimen data. Below is an overview of studies by topic.

* private studies, ones which are uploaded and marked private, are not counted here

News

October 02, 2014: ICBI Symposium 2014

[read]


May 02, 2014: Featured in Frontiers' Top 10 2013 Most viewed Genetics Research articles [read]


March 12, 2014: AAAS Big Data Blog [read]





G-DOC Plus Launch Pad!


Welcome! The G-DOC Plus Launch Pad is your one-stop resource for learning more about G-DOC and getting started on the platform.











G-DOC has over seventy studies We know this can be overwhelming! Let us guide you to choose the study that is relevant for your research.

Let's Go! >

Help


kb472 ¬

What's your area of interest?

G-DOC Plus has three overlapping entry points for the user based on their interests. Choose your area of interest to launch the workflow.















Select disease/data of interest

DEMENTIA

1 30 30 study samples biospecimen

PEDIATRIC CANCERS

5 256 256 biospecimen

BREAST CANCER

25 3952 4532 studies samples biospecimen

COLON CANCER

10 1226 1262 studies samples biospecimen

LUNG CANCER

1 478 443 study samples biospecimen

MUSCULAR DYSTROPHY

1 36 36 biospecimen

LIVER CANCER

3 298 468 studies samples biospecimen

CELL_LINE_COLLECTIONS

1 60 59 study samples biospecimen

OVARIAN CANCER

1 1711 564 biospecimen

PROSTATE CANCER

1 465 538 study samples biospecimen

PANCREATIC CANCER


1 52 51 study samples biospecimen

BRAIN CANCER


1 671 804 samples biospecimen



Choose between patient and cell line data




Select study




Note: If you click on "More", you will see complete description of the dataset


Study selected. Let's now create some groups



Explore clinical data or create groups



Let's save the Oligodendroglioma patients with tumor grade II and grade III



Clinical evaluation

Karnofsky @


Disease evaluation by MRI @

Neurologic exam score @ Performance Status Score:

Repeat this procedure to create as many number of patient groups as needed.

Note: "View detailed report" shows a detailed clinical report of those selected patients.


Enter a name for the list, and click "Save"



Perform Group Comparison, then perform Gene expression KM analysis

METHOD 1

Now we want to compare the Oligodendroglioma tumor grade II patients with Oligo grade III patients





Perform Group Comparison Analysis

Current Study: REMBRANDT

Select a baseline group and a comparison group(s)

Rem_olig_grade2	\$.05	False Discovery Rate(FDR): Bt 💠
Select comparison group		Fold Change	Data-Type
Rem_olig_grade3	•	1.5	GENE EXPRESSION \$
		Statistical Method	Dataset
		T-Test: Two Sample Test \$	mas5 normalization \$

Submit Analysis

Select baseline group (less screwed up group), comparison group, and settings for the comparison analysis.

Click "Submit analysis"



Home

Studies

Lists

Analyse

Groups

Notifications

Study Options -

Help



kb472 ~

Notifications

Below are your latest running analyses. Once completed, click on the Analysis name to see detailed results.

GROUP_COMPARISON (2:05 1/22/2015)

Complete



Once you click "submit analysis", you will be re-directed to the "Notifications" page. Once the status is "Complete", click on "Group comparison"

Results of group comparison



Home

Studios

Lists

Analyse

Groups

Notifications

Study Options


Help

Q

kb472 ~

Analysis Results

Current Study: REMBRANDT



List Name:


Save Selected ↓

View HeatMap for selected reporters


These are the results of the group comparison.

You can sort this table based on any of the columns.


In this example, they are sorted based on fold change



Going from Group Comparison to Gene expression KM



Gene expression KM plot




Create group of interest, perform gene expression KM

METHOD 2

Steps

- Login
- Select study of interest (follow steps in slides 4–8)
- Create group of interest (follow steps in slides 9-12)
- Go to Study Options Select "Gene expression KM plot"
- Select settings based on interest, click "plot"

See screenshots on next slide





Current Study: REMBRANDT change study?

Select Patient Groups: Rem_oligo_gr2

Select Gene: PAPPA


Select Endpoint: Overall survival in months

Data-Type: GENE EXPRESSION

Dataset: mas5 normalization

Plot Reset

***....performing gene expression analysis





Perform gene expression KM plot on all samples


METHOD 3

Steps

- Login
- Select study of interest (follow steps in slides 4–8)
- Go to "Study Options" Select "Gene expression KM plot"
- Select settings based on interest, click "plot"
 See screenshots on next slide







Changing settings for your KM plot

- No matter which method you use to create your KM plots, you can
 - Change "fold change" values and "Redraw" your plots to see how the patient stratification changes
 - Choose another "reporter id" (gene expression probe id) and "Redraw" your plots to see how the patient stratification changes
 - Save your plot click on "Chart options" and you should

Algorithm used in gene expression KM plot tool

- There are various ways of defining "low", "intermediate", and "high" expression. Some methods include: using mean or median as cut off; using 25% and 75% quartiles, etc.
- In G-DOC, we use "fold change" concept to split the gene expression data into "low", "intermediate" and "high expression".
 - The fold change is calculated for each sample in the input based on mean expression of a gene (indicating how many fold higher or lower the expression is compared to mean).
 - G-DOC lets users to choose the fold change cut off values. So if the cut off is 2, then samples with fold change < -2 are equivalent to "low expression" (down-regulated). Samples with fold change between -2 and 2 are equivalent to "intermediate expression". Samples with fold change > 2 are equivalent to "high expression" (up-regulated).
- Note that the gene expression values are log2 values, so fold change 1 means ratio of 2.

TIPS

General tips

 G-DOC Plus works best is you don't use the back button in the web browser repeatedly.

Once you select a study, most tools will be easily available from the top menu bar inside G-DOC *Plus*.

 The <u>Pathway enrichment</u> and the <u>Lists</u> tool may sometimes take a few seconds longer to execute than other tools (since they are directly connecting to the server every time). Your patience is highly appreciated.

Clearing cache

- If the G-DOC web page does not respond after several seconds, try:
 - refreshing the page.
 - Log out and log back in, and try again
 - If the above two do not work, its possible that your web browser cache may need to be cleared
 - For Google chrome, go to Settings -> Show Advanced Settings -> Under "Privacy", select Clear Browsing data

For Mozilla Firefox, go to Preferences -> Advanced -> Network -> Under "Cached Web

Content" -> Clear now





We are working hard to improve G-DOC Plus.
 Please feel free to email your questions and comments (no homework questions please) to us at: gdoc-help@georgetown.edu