

How to grid enable prostate SPORES DE caTissue Suite v1.1.2 - draft



To Print the Guide

You can create a PDF of the guide. For instructions refer to the tip [How do I print multiple pages?](#) If you want to print a single page, refer to [How do I print a page?](#)

The *How to grid enable prostate SPORES DE* has the following sections.

- [Copyright Notice](#)
- [License](#)
- [About This Guide](#)
 - [Purpose](#)
 - [Audience](#)
 - [Related Documentation](#)
 - [Document Change History](#)
 - [Contacts and Support](#)
- [Importing the DE form](#)
- [Exporting the DE form](#)
- [Generating the caCORE API service](#)
- [Running Introduce to create the caGRID service](#)
- [Running the Dynamic Extensions caGrid Test Queries](#)

Copyright Notice

Copyright 2010 Washington University in St. Louis ("caBIG® Participant"). caTissue was created with NCI funding and is part of the caBIG® initiative. The software subject to this notice and license includes both human readable source code form and machine readable, binary, object code form (the "caBIG® Software").

This caBIG® Software License (the "License") is between caBIG® Participant and You. "You (or "Your") shall mean a person or an entity, and all other entities that control, are controlled by, or are under common control with the entity. "Control" for purposes of this definition means (i) the direct or indirect power to cause the direction or management of such entity, whether by contract or otherwise, or (ii) ownership of fifty percent (50%) or more of the outstanding shares, or (iii) beneficial ownership of such entity.

License

Provided that You agree to the conditions described below, caBIG® Participant grants You a non-exclusive, worldwide, perpetual, fully-paid-up, no-charge, irrevocable, transferable and royalty-free right and license in its rights in the caBIG® Software, including any copyright or patent rights therein, to (i) use, install, disclose, access, operate, execute, reproduce, copy, modify, translate, market, publicly display, publicly perform, and prepare derivative works of the caBIG® Software in any manner and for any purpose, and to have or permit others to do so; (ii) make, have made, use, practice, sell, and offer for sale, import, and/or otherwise dispose of caBIG® Software (or portions thereof); (iii) distribute and have distributed to and by third parties the caBIG® Software and any modifications and derivative works thereof; and (iv) sublicense the foregoing rights set out in (i), (ii), and (iii) to third parties, including the right to license such rights to further third parties. For sake of clarity, and not by way of limitation, caBIG® Participant shall have no right of accounting or right of payment from You or Your sublicensees for the rights granted under this License. This License is granted at no charge to You. Your downloading, copying, modifying, displaying, distributing or use of caBIG® Software constitutes acceptance of all of the terms and conditions of this Agreement. If you do not agree to such terms and conditions, you have no right to download, copy, modify, display, distribute or use the caBIG® Software.

1. Your redistributions of the source code for the caBIG® Software must retain the above copyright notice, this list of conditions and the disclaimer and limitation of liability of Article 6 below. Your redistributions in object code form must reproduce the above copyright notice, this list of conditions and the disclaimer of Article 6 in the documentation and/or other materials provided with the distribution, if any.
2. Your end-user documentation included with the redistribution, if any, must include the following acknowledgment: "This product includes software developed by Washington University in St. Louis." If You do not include such end-user documentation, You shall include this acknowledgment in the caBIG® Software itself, wherever such third-party acknowledgments normally appear.
3. You may not use the names "Washington University in St. Louis", "The National Cancer Institute", "NCI", "Cancer Bioinformatics Grid" or "caBIG®" to endorse or promote products derived from this caBIG® Software. This License does not authorize You to use any trademarks, service marks, trade names, logos or product names of either caBIG® Participant, NCI or caBIG®, except as required to comply with the terms of this License.
4. For sake of clarity, and not by way of limitation, You may incorporate this caBIG® Software into Your proprietary programs and into any third party proprietary programs. However, if You incorporate the caBIG® Software into third party proprietary programs, You agree that You are solely responsible for obtaining any permission from such third parties required to incorporate the caBIG® Software into such third party proprietary programs and for informing Your sublicensees, including without limitation Your end-users, of their obligation to secure any required permissions from such third parties before incorporating the caBIG® Software into such third party proprietary software programs. In the event that You fail to obtain such permissions, You agree to indemnify caBIG® Participant for any claims

against caBIG® Participant by such third parties, except to the extent prohibited by law, resulting from Your failure to obtain such permissions.

- For sake of clarity, and not by way of limitation, You may add Your own copyright statement to Your modifications and to the derivative works, and You may provide additional or different license terms and conditions in Your sublicenses of modifications of the caBIG® Software, or any derivative works of the caBIG® Software as a whole, provided Your use, reproduction, and distribution of the Work otherwise complies with the conditions stated in this License.
- THIS caBIG® SOFTWARE IS PROVIDED "AS IS" AND ANY EXPRESSED OR IMPLIED WARRANTIES (INCLUDING, BUT NOT LIMITED TO, THE IMPLIED WARRANTIES OF MERCHANTABILITY, NON-INFRINGEMENT AND FITNESS FOR A PARTICULAR PURPOSE) ARE DISCLAIMED. IN NO EVENT SHALL WASHINGTON UNIVERSITY IN ST. LOUIS OR ITS AFFILIATES BE LIABLE FOR ANY DIRECT, INDIRECT, INCIDENTAL, SPECIAL, EXEMPLARY, OR CONSEQUENTIAL DAMAGES (INCLUDING, BUT NOT LIMITED TO, PROCUREMENT OF SUBSTITUTE GOODS OR SERVICES; LOSS OF USE, DATA, OR PROFITS; OR BUSINESS INTERRUPTION) HOWEVER CAUSED AND ON ANY THEORY OF LIABILITY, WHETHER IN CONTRACT, STRICT LIABILITY, OR TORT (INCLUDING NEGLIGENCE OR OTHERWISE) ARISING IN ANY WAY OUT OF THE USE OF THIS caBIG® SOFTWARE, EVEN IF ADVISED OF THE POSSIBILITY OF SUCH DAMAGE.

About This Guide

This document details a specific scenario where a custom dynamic extension form is grid enabled for caTissue Suite 1.1.2.

Purpose

The *How to grid enable prostate SPORES DE* provides details for the following:

- Importing the DE form
- Exporting the DE form
- Generating the caCORE API service
- Correcting the xmi file with Enterprise Architect
- Running Introduce to create the caGRID data service
- Deploying the Dynamic Extension grid service
- Running test queries

Audience

This guide is intended for technical IT staff. It describes steps involved in importing dynamic extension forms into caTissue and deploying the dynamic extension caGrid service which involves deploying multiple JBoss containers on the same machine, as well as configuring multiple system settings.

Related Documentation

You can find the latest documentation for the caTissue Suite at the [caTissue Documentation Wiki](#)

Document Change History

Version Number	Date	Description	Contributor
1.0	7/25/2011	New document	TBTT Knowledge Center

Contacts and Support

Knowledge Center Main Page	https://cabig-kc.nci.nih.gov/Biospecimen/KC/index.php/Main_Page
Knowledge Center Discussion Forum	https://cabig-kc.nci.nih.gov/Biospecimen/forums/
Knowledge Center Contact	tbpt_kc_support@mga.wustl.edu

Importing the DE form

The instructions given below are abbreviated. For more details on the import process please refer to [Chapter 5 of the caTissue Suite 1.1.1 Technical Manual](#)

- Open a command prompt and change the directory to where the caTissue package was unzipped
- From this directory run the following ant target: `ant -f deploy.xml import_xmi -Dfilename="<path>VPSBINParticipant-3.xmi" -DhookEntity="<HookEntityClassName>" -DmainContainerList="<MainContainerCSVFileName>" -Dpackage="<Package>" -Dcondition="<ConditionCSVFileName>"`
- Note that the correct path needs to be set along with the proper separator "/" or "\".

Exporting the DE form

Video demonstration of the following steps is available [here](#)

The instructions given below are abbreviated. For more details on the import process please refer to [Chapter 5 of the caTissue Suite 1.1.1 Technical Manual](#)

1. Open a command prompt and change the directory to where the caTissue package was unzipped
2. From this directory run the following ant target: `ant -f deploy.xml export_xmi -Dgroupname=PSBINParticipant-3 -Dfilename=<your_path>\PSBINParticipant-3.xmi -Dversion 1.1`

Generating the caCORE API service

These instructions are taken from a [FAQ](#) posted on the TBTT Knowledge Center

1. Copy and unzip the caCORE SDK 3.2 from the following location http://gforge.nci.nih.gov/frs/download.php/1938/caCORE_SDK_321.zip
2. Copy the output file from the *Exporting the DE Form* step to the caCORE_SDK_HOME\models\xmi directory.
3. Modify the following properties in caCORE_SDK_HOME\conf\deploy.properties:
 - a. System Properties
 - i. j2se_container_home should point to your local install
 - b. Project Properties
 - i. project_name=PSBINParticipant
 - ii. webservice_name=PSBINParticipant
 - c. Model Properties
 - i. model_filename=PSBINParticipant-3.xmi
 - ii. fixed_filename=fixed_PSBINParticipant-3.xmi
 - iii. include_package=.*ClinicalAnnotation.*
 - d. Database Properties (e.g. db properties as set in the caTissueInstall.properties file)
 - e. J2SE CONTAINER PROPERTIES, Set "start_container=no"
 - f. Set "disable_writable_api_generation=no".
 - g. For more information on caCORE, refer to the caCORE v3.2 User Guide ftp://ftp1.nci.nih.gov/pub/cacore/SDK/v3.2.1/caCORE_SDK_3.2.1_Programmers_Guide.pdf
4. Navigate to the caCORE_SDK_HOME directory
5. Run "ant build-system"
6. Modify the following two properties from caTissueInstall.properties
 - i) caCORE.jBoss.home.dir as path of the jboss where caTissue is deployed.
 - ii) [caCORE.project.name](#) should be same as project_name you specified deploy.properties
7. Go to command prompt and navigate to <caTissue_installable>\catissue_de_integration_client
8. Run ant copyDeIntegrationWar.
9. Run ant modifyProjectWar.

Running Introduce to create the caGRID service

The caGRID 1.2 software is a prerequisite for this step. You can download the software from this [link](#) and view installation documentation [here](#).

1. Create a new directory for the caGRID service (the output from step 3.b will be saved here)
2. From the caGrid installation directory run `ant introduce`
3. **See *Correcting the xmi file with Enterprise Architect before proceeding***
4. Using the Introduce: Grid Service Authoring Toolkit
 - a. Select the Create caGrid Service Skeleton from the menu bar
 - b. Introduce STEP 1 - Browse to the directory created in step #1 above
 - c. Introduce STEP 2 - Enter the name of the service as "PSBINParticipant"
 - d. Introduce STEP 3 - Enter a Java Package for the generated code as "PSBINParticipant"
 - e. Introduce STEP 4 - No action
 - f. Change the Service from Analytical to Data
 - g. Click on the Create button
 - h. From the Data Service Configuration dialog box, select the caCORE SDK v 3.2.(1) option and click OK
 - i. From the caCORE SDK v 3.2.(1) Style: Create caCORE SDK Backend caGRID Data Service wizard click Next: Client Selection
 - j. Click Browse to select the Client Lib Directory. This is located in the <your_caCORE_install>\output\PSBINParticipant\package\client\lib directory
 - k. Click Next: Configuration
 - l. In the Remote Service URL enter the URL for caTissue. For example: <http://localhost:8080/catissuecore>
 - m. Click Next: Model
 - n. **See *Correcting the xmi file with Enterprise Architect before proceeding to the next step***
 - o. Select Domain Model from File and navigate to the file that was output from the *Correcting the xmi file with Enterprise Architect* step using the Browse button
 - p. Select "Fix EA Model" check box and then click on Browse button and browses to the <your_caCORE_install>.
 - q. Enter a Project Short Name: PSBINParticipant
 - r. Enter a Project Version: 1.0
 - s. Click OK
 - t. You might see an Errors dialog at this point related to GME, caDSR and errors parsing the XMI file. It is acceptable to ignore the

- error by clicking on the Hide button
- u. A pop up will appear. Click on Resolve, a pop up will appear select the xsd from
<your_caCORE_install>\output{color:#000000}PSBINParticipant\package\client\conf. Click on ok.
Click Next: Schemas
Click the Done button
The Modify Service Interface dialog will be displayed
Click the Save button and Confirm the save.
This completes the service creation.

h3. Correcting the xmi file with Enterprise Architect

If this step is not done then running introduce will throw an error

- # Start Enterprise Architect
- # Create a new project called PSBINParticipant.eap and click the Save button
- # From the Select Model(s) dialog check Domain Model and click the OK button
- # Select Project from the main menu bar and go to Import/Export -> Import Package from XMI
- # Navigate to the XMI file exported from caTissue in the Exporting the DE Form section. (e.g. <your_path>\PSBINParticipant-3.xmi) and click the Import button
- # Select Yes if prompted to overwrite the current package
- # Select Close when the import is complete
- # From the Project Browser view right-click on the Logical View
- # Navigate to Import/Export -> Export package to XMI file
- # From the Export Package to XMI dialog select the desired output directory, enter a value for the File name and click Save.
- # From the Export Package to XMI dialog click on Export and then Close.

h3. Deploying Dynamic Extension Grid Services

caTissue Suite v1.2 is grid-enabled wherein all the caGrid queries are routed through a single designated caTissue user account, for example, a caGrid user will not be directly provisioned in caTissue. In a future release of caTissue, one will be able to provision the caGrid users in caTissue.



Caution

The single designated caTissue user should be a user with the role Scientist who does not have access to any identified data AND is not a Principal Investigator or Coordinator of any collection protocol. The caTissue instance to which the grid service is going to point, should be deployed on a HTTP JBoss instance. In this release, the grid service is not supported if caTissue is deployed on a secure JBoss instance. Please refer to [Chapter 5 - Private Public Data Store](#) of the *caTissue Suite Deployment Guide* for details on how to de-identify caTissue data.

h4. Background

The caTissue caGrid data service requires that a caGrid service (or grid node) is deployed and functioning. Below are detailed instructions for deploying the service but assumptions are made that might not fit your environment or deployment needs. For more information or assistance with deploying the caGrid service, please see the [caGrid Knowledge Center](#). Multiple data services can use one grid node but require separate JBoss containers. Please see the [TBTT FAQ](#) on changing the JBoss default ports to run multiple JBoss containers from the same machine. By changing the default ports you can deploy multiple versions of JBoss without port conflicts.

The deployment steps for the caGrid Service are the following:

- # Deploy JBoss 4.0.4.
- # Deploy Java JDK 1.5.
- # Deploy the caGrid 1.2 service (using the caGrid installer will automatically download Ant and Globus).
- # Create the secure caGrid service.

The deployment steps for the caTissue data services are the following:

- # Deploy Globus into the JBoss for the desired caTissue service.
- # Deploy and configure the desired caTissue service.
- # Start the caTissue service JBoss.

h4. Prerequisites

**Note**

This section is for reference only. This guide assumes a functioning caGrid service deployment on the same machine where one of the caTissue dynamic extension grid services will be deployed. See the [caTissue Suite 1.2 Deployment Guide - Chapter 6, Deploying caTissue caGRID Data Service v1.2](#) for more details about deploying a caGrid service.

```

|| Software Name || Version || URL || Install Guide ||
| caGrid | 1.2 | http://gforge.nci.nih.gov/frs/download.php/3738/caGrid-installer-1.2.zip |
http://gforge.nci.nih.gov/plugins/scm cvs/cvsweb.php/checkout
/cagrid-1-0/Documentation/docs/installer/caGrid-1-2_Installer_Guide.pdf?rev=HEAD;content-type=application%2Foctet-stream;cv
|
| Globus Toolkit | 4.0.3 | http://gforge.nci.nih.gov/frs/download.php/1334/ws-core-enum-4.0.3.zip |
http://www.globusconsortium.org/tutorial/ch6/ |
| Jboss | 4.0.4 GA | http://labs.jboss.com/jbossas/downloads/ | |

```

h4. Deploying and Securing Globus

Follow the directions [here](#) for securing Globus.

h4. Deploying the Data Service

**Note**

Dynamic Extensions model data services should be deployed on a different grid JBoss than the caTissue grid JBoss instance. If the caTissue service & Dynamic Extensions services are deployed on the same grid JBoss, then only one of the services will work.

```

# Files that were created from the Running Introduce to create the caGRID service step will be referenced as
introduce_output. Copy the credentials.properties file located in the
CATISSUE_HOME/CaTissueSuite_caGrid_Server to the introduce_output folder.
# Open the command prompt, and change the directory to the introduce_output.
# Update the file credentials.properties. (copied from step 1)
|| Parameter Name || Details ||
| user.name | caTissue user name to which the caGrid users should be mapped.

```

**Caution**

The Administrator should ensure that the configured user configured should have a Scientist role in the application. That user should not be Principal Investigator or Coordinator of any collection protocol.

**Note**

This caTissue account will expire as per the property settings in `password.expire_after_n_days`. It is highly recommended that you set a calendar alert to notify administrators that this account password will expire in *n* days. The grid service will not be accessible if the account expires.

|
| `password` | Password of the common user to which all the caGrid users will be mapped. |
| `keystoreFilePath` | Description: Path to the keystore file if the caTissue is deployed as HTTPS. This key should match to the key of caTissue web application instance. The path must be separated by / (forward slash) character. To generate the keystore file, perform the first step mentioned in the section [Configuring JBoss Server to deploy caTissue using HTTPS](#).
| **Default Value:** N/A
| **Permissible Value:** N/A |

- After editing the file, copy it to the following: `{user.home}/catissueservice` where `{user.home}` ...which is the Windows default user directory, for example, `C:\Documents and Settings\srikanth\ catissueservice`. In Unix or Linux, it is the user home directory, for example, `/home/srikanth/catissueservice`.

**Note**

This file will have to be copied to the user home of the user under which the JBoss server is going to be run.

- Update the file `serviceMetadata.xml` as follows:
 - Location:** `introduce_output/etc`
 - Locate ns3:** `PointOfContact` affiliation and add in the details for email, firstName, lastName, and role. `PhoneNumber` is optional.
 - Locate ns9:** `ResearchCenter`, `displayName`, and update details such as your research center `displayName`, `shortName`, and address information. These details will be used to place the service on the Google map on the [Grid Portal Page](#). Details about the person maintaining the service should also be added under `PointOfContact`.
- Update the file `service.properties` as follows:
 - ** *Location:** `introduce_output/`
The property `cqlQueryProcessorConfig_appserviceUrl` should point to the host and port of the caTissue web application instance in the following format: `http://host:port/CA/http/remoteService`
...where the `host` and `port` are the hostname and port number of the caTissue web application, for example, `http://testserver3.wustl.edu:8080/PSBINParticipant/http/remoteService`.
- Update the file `introduce_output/build-deploy.xml`.
 - Search for line: `<property name="jboss.dir" value="${env.JBOSS_HOME}"/>`
 - Replace `${env.JBOSS_HOME}`
...with the folder name of `JBOSS_HOME`, for example, `<property name="jboss.dir" value="/usr/local/jboss-4.0.4.GA/"`
 - From this folder, run the following commands: `>ant clean`
`>ant all`
`>ant deployJBoss`
- Start the JBoss server.
 - Open a command prompt.
 - Change the directory `cd %JBOSS_HOME%\bin`.
 - Run the command `run.bat -c default`.

**Note**

You can verify that the grid service has been successfully deployed by typing the following URL in the browser:
`https://<hostname>:<port>/wsrf/services/cagrid/PSBINParticipant`

5. To stop the JBoss server:
 - a. Open a command prompt.
 - b. Change the directory `cd %JBOSS_HOME%\bin.`
 - c. Run the command `shutdown.bat -S.`
Similarly other Dynamic Extensions & Clinical Annotations model CaGrid services can be deployed by following the same steps from corresponding folder.

Running the Dynamic Extensions caGrid Test Queries

1. Log in to the grid portal (where the grid node was pointed)
 - a. [Training portal](#)
 - b. [Production portal](#)
2. Go to the [query builder page](#)
3. Insert your grid service URL (under the select UML class tab)
- 4.