

# caNanoLab 2.3.6 Release Note

Welcome to the caNanoLab 2.3.6 Release Note, January 7, 2019. Topics in this Release Note include the following.

- [Improvements](#)
- [Known Issues and Workarounds](#)
- [Release History](#)



## JIRA

If you have access to the JIRA Tracking System, click a JIRA link to review additional development details.

## Improvements

The following table lists and describes the new feature(s) in this caNanoLab release.

Description of the improvement(s)	JIRA issue link
Added the % dose per g tissue unit for pharmacokinetics, biodistribution, and tumor accumulation assays.	<a href="#">CANANO LAB-568</a>
For composition, when you save/update a composing entity, caNanoLab now consistently displays the amount(s) saved. Previously, you may only have seen amount(s) when editing the entity but not after saving/updating.	<a href="#">CANANO LAB-571</a>

## Known Issues and Workarounds

The following table lists and describes how to work around the following remaining issues in caNanoLab.

Known Issue	Workaround	JIRA issue link
You cannot add multiple points of contact to one organization.	Add a division or department name to the organization's name and then add the contact to that organization.	<a href="#">CANAN OLAB-61</a>
The Results tab is populated when a curator/administrator submits a batch job to create a data availability matrix, but the tab disappears if the curator/administrator logs out and logs in again, even if the job is still running.	None	<a href="#">CANAN OLAB-81</a>
When submitting a new publication, clicking Reset does not reset the Submit Publication Form.	None	<a href="#">CANAN OLAB-197</a>
You cannot sort Protocol search results by the Protocol Type and Protocol Name columns.	None	<a href="#">CANAN OLAB-367</a>
You have to click the DOI external URL twice (or open the link in a new window by right-clicking) to view the external PubMed content.	None	<a href="#">CANAN OLAB-368</a>
A Publication Search returns more records than is valid based on the search criteria.	None	<a href="#">CANAN OLAB-391</a>

## Release History

Release	Date
2.3.6	7 January 2019
2.3.5.1	15 November 2018
2.3.5	2 October 2018

2.3.4	13 July 2018
2.3.3	25 January 2018
2.3.2	October 2017
2.3	April 2017
2.2	January 2017
2.1.1	February 2016
2.1	November 2015
2.0	December 2014
1.5.5	May 2013
1.5.4	August 2012
1.5.3	May 2012
1.5.2	October 2010