

CTP Properties

Life Sciences Distribution (LSD) 1.2 Installation Wizard: ncia

Installation Options

Step 32 of 46

CTP Server Configuration

National Cancer Institute

WARNING WARNING WARNING WARNING WARNING WARNING WARNING WARNING

The installer (for install or upgrade) will remove any previously installed application servers (using caBIG installers) based on the information you provide.

If you are doing an upgrade and the CTP server stores images in the default relative location (storage), they could be lost. Please be sure to copy your image files out of this relative location before doing the upgrade.

If you are doing a fresh install, it is recommended you select an absolute location for storage that is NOT under the NBIA installation directory.

CTP Server HTTP port: 23080

CTP HTTP Import Service port: 23104

CTP Server Storage Service port: 23204

CTP Verifier Service port: 23304

CTP Server Admin Password:

Confirm Password:

CTP Server King Password:

Confirm Password:

Directory for Storage Service: storage

Index directory for Storage Service: index

☐ Require authentication on Import Service

Copyright (c) 2008 National Cancer Institute - <http://www.cancer.gov>

Previous Next Quit

- ctp.server.port=23080
- ctp.server.http_import_service.port=23104
- ctp.server.storage_service.port=23204
- ctp.server.verifier.port=23304
- ctp.server.admin.password=password2
- ctp.server.king.password=password2
- ctp.server.storage_service.root=storage
- ctp.server.storage_service.index=index

Life Sciences Distribution (LSD) 1.2 Installation Wizard: ncia

Installation Options

Step 33 of 46



CTP Client Configuration

CTP Client HTTP port:	24080
CTP Client HTTP Import Service port:	24104
CTP Client DICOM Import Service port:	24304
CTP Client Export Service port:	24204
CTP Client Admin Password:	*
Confirm Password:	*
CTP Client King Password:	*
Confirm Password:	*

Copyright (c) 2008 National Cancer Institute - <http://www.cancer.gov>

Previous Next Quit

- ctp.server.require_authentication=no
- ctp.client.port=24080
- ctp.client.http_import_service.port=24104
- ctp.client.dicom_import_service.port=24304
- ctp.client.export_service.port=24204
- ctp.client.admin.password=password2
- ctp.client.king.password=password2