

# caNanoLab 2.3.7 Release Note

Welcome to the caNanoLab 2.3.7 Release Note, June 11, 2019. Topics in this Release Note include the following.

- [Improvements](#)
- [Fixed Bugs](#)
- [Known Issues and Workarounds](#)
- [Release History](#)





## JIRA

If you have access to the JIRA Tracking System, click a JIRA link to review additional development details.



## Improvements

The following table lists and describes the improvement(s) in this caNanoLab release.

| Description of the improvement(s)   | JIRA issue link               |
|---|-------------------------------|
| <p>To simplify the process of adding data to caNanoLab, you can now import tables of data from, for example, Microsoft Excel, into the Findings section of assay data entry.</p> <ol style="list-style-type: none"><li>1. Edit a Sample.</li><li>2. Select the Characterization menu, and add/edit a characterization.</li></ol>  <ol style="list-style-type: none"><li>3. Add/Edit a Finding, and click <b>Import csv</b> to add data from a file.</li></ol>  <ol style="list-style-type: none"><li>4. Click <b>Save</b>.</li></ol> | <a href="#">CANANOLAB-570</a> |

## Fixed Bugs

The following table lists and describes the bug(s) resolved in this caNanoLab release.

| Description of the bug fix(es)  | JIRA issue link               |
|---|-------------------------------|
| <p>If you add an Inherent Function in a Composing Element, the Inherent Function is now saved.</p> <ol style="list-style-type: none"><li>1. Edit a Sample.</li><li>2. Select the Composing Element menu and add an Inherent Function.</li></ol>   <ol style="list-style-type: none"><li>3. Click <b>Save</b> and the Inherent Function is saved in the Composing Element.</li></ol> | <a href="#">CANANOLAB-481</a> |
| <p>When you edit or add a characterization (CharacterizationFindingFinding InfoData and Conditions), changing the order of the columns does not delete all the data.</p>  | <a href="#">CANANOLAB-584</a> |
| <p>When you save Finding Info data and there is a server error, you will receive an error message indicating that your data was not saved.</p>  | <a href="#">CANANOLAB-585</a> |

|   |                               |
|---|-------------------------------|
| When you edit or add a characterization (CharacterizationFindingFinding InfoData and Conditions), the text changes to red if your entry causes an error.                                | <a href="#">CANANOLAB-586</a> |
| The correct help page now appears on the Publication page, <a href="#">Managing Publications</a> .  | <a href="#">CANANOLAB-594</a> |
| There is now one help button on the Create Publication page, <a href="#">Create a Publication</a> .   | <a href="#">CANANOLAB-595</a> |
| If you log into caNanoLab as a Curator and generate data availability for all samples, the correct help page now appears on the Results page, <a href="#">Reviewing Batch Results</a> . | <a href="#">CANANOLAB-596</a> |

## Known Issues and Workarounds

The following table lists and describes how to work around the following remaining issues in caNanoLab.

| Known Issue   | Workaround  | JIRA issue link               |
|---|---|-------------------------------|
| You cannot add multiple points of contact to one organization.  | Add a division or department name to the organization's name and then add the contact to that organization. | <a href="#">CANANOLAB-61</a>  |
| The Results tab is populated when a curator/administrator submits a batch job to create a data availability matrix, but the tab disappears if the curator/administrator logs out and logs in again, even if the job is still running. | None  | <a href="#">CANANOLAB-81</a>  |
| When submitting a new publication, clicking Reset does not reset the Submit Publication Form.   | None  | <a href="#">CANANOLAB-197</a> |
| You cannot sort Protocol search results by the Protocol Type and Protocol Name columns.   | None  | <a href="#">CANANOLAB-367</a> |
| You have to click the DOI external URL twice (or open the link in a new window by right-clicking) to view the external PubMed content.  | None  | <a href="#">CANANOLAB-368</a> |
| A Publication Search returns more records than is valid based on the search criteria.   | None  | <a href="#">CANANOLAB-391</a> |

## Release History

| Release | Date             |
|---------|------------------|
| 2.3.7   | 11 June 2019     |
| 2.3.6   | 7 January 2019   |
| 2.3.5.1 | 15 November 2018 |
| 2.3.5   | 2 October 2018   |
| 2.3.4   | 13 July 2018     |
| 2.3.3   | 25 January 2018  |
| 2.3.2   | October 2017     |
| 2.3     | April 2017       |
| 2.2     | January 2017     |
| 2.1.1   | February 2016    |
| 2.1     | November 2015    |
| 2.0     | December 2014    |
| 1.5.5   | May 2013         |
| 1.5.4   | August 2012      |

|       |              |
|-------|--------------|
| 1.5.3 | May 2012     |
| 1.5.2 | October 2010 |