

caNanoLab 2.3.7 Release Note

Welcome to the caNanoLab 2.3.7 Release Note, June 11, 2019. Topics in this Release Note include the following.

- [Improvements](#)
- [Fixed Bugs](#)
- [Known Issues and Workarounds](#)
- [Release History](#)

 **JIRA**
If you have access to the JIRA Tracking System, click a JIRA link to review additional development details.

Improvements

The following table lists and describes the improvement(s) in this caNanoLab release.

Description of the improvement(s)	JIRA issue link
To simplify the process of adding data to caNanoLab, you can now import tables of data from, for example, Microsoft Excel, into the Findings section of assay data entry. 1. Edit a Sample. 2. Select the Characterization menu, and add/edit a characterization. GENERAL INFO COMPOSITION CHARACTERIZATION PUBLICATION 3. Add/Edit a Finding, and click Import csv to add data from a file. Finding Add and Import csv 4. Click Save.	CANANOLAB-570

Fixed Bugs

The following table lists and describes the bug(s) resolved in this caNanoLab release.

Description of the bug fix(es)	JIRA issue link
If you add an Inherent Function in a Composing Element, the Inherent Function is now saved. 1. Edit a Sample. 2. Select the Composing Element menu and add an Inherent Function. GENERAL INFO COMPOSITION CHARACTERIZATION PUBLICATION and Inherent Function Add 3. Click Save and the Inherent Function is saved in the Composing Element.	CANANOLAB-481
When you edit or add a characterization (CharacterizationFindingFinding InfoData and Conditions), changing the order of the columns does not delete all the data.	CANANOLAB-584
When you save Finding Info data and there is a server error, you will receive an error message indicating that your data was not saved.	CANANOLAB-585

When you edit or add a characterization (CharacterizationFindingFinding InfoData and Conditions), the text changes to red if your entry causes an error.	CANANOLAB-586
The correct help page now appears on the Publication page, Managing Publications .	CANANOLAB-594
There is now one help button on the Create Publication page, Create a Publication .	CANANOLAB-595
If you log into caNanoLab as a Curator and generate data availability for all samples, the correct help page now appears on the Results page, Reviewing Batch Results .	CANANOLAB-596

Known Issues and Workarounds

The following table lists and describes how to work around the following remaining issues in caNanoLab.

Known Issue	Workaround	JIRA issue link
You cannot add multiple points of contact to one organization.	Add a division or department name to the organization's name and then add the contact to that organization.	CANANOLAB-61
The Results tab is populated when a curator/administrator submits a batch job to create a data availability matrix, but the tab disappears if the curator/administrator logs out and logs in again, even if the job is still running.	None	CANANOLAB-81
When submitting a new publication, clicking Reset does not reset the Submit Publication Form.	None	CANANOLAB-197
You cannot sort Protocol search results by the Protocol Type and Protocol Name columns.	None	CANANOLAB-367
You have to click the DOI external URL twice (or open the link in a new window by right-clicking) to view the external PubMed content.	None	CANANOLAB-368
A Publication Search returns more records than is valid based on the search criteria.	None	CANANOLAB-391

Release History

Release	Date
2.3.7	11 June 2019
2.3.6	7 January 2019
2.3.5.1	15 November 2018
2.3.5	2 October 2018
2.3.4	13 July 2018
2.3.3	25 January 2018
2.3.2	October 2017
2.3	April 2017
2.2	January 2017
2.1.1	February 2016
2.1	November 2015
2.0	December 2014
1.5.5	May 2013
1.5.4	August 2012

1.5.3	May 2012
1.5.2	October 2010