

calnt 0005 - Adding a new array platform in calIntegrator

Question: How do I add a new array platform in calIntegrator?

Topic: calIntegrator Usage

Release: calIntegrator 1.2 and above

Date entered: 03/23/2011

Answer

You can create a study without genomic data, but you cannot add genomic data to a calIntegrator study without a corresponding supported array platform. Currently, calIntegrator only supports following types of platforms: Affymetrix Gene Expression, Affymetrix Copy Number, Agilent Gene Expression, and Agilent Copy Number. Please refer to the [User Guide](#) for step-by-step instructions on how to create a new platform.

To add a new platform, the array annotation file is required. For Affymetrix platforms, please refer to [GeWorkbench025](#) for how to obtain the annotation file. For other platforms, please refer to the following resources:

- [caArray](#)
- [GEO](#)
- [Arrayexpress](#)

Have a comment?

Please leave your comment in the [calIntegrator End User Forum](#).