

calnt 0008 - An import of an Agilent array dataset from caArray into calIntegrator 2 failed

Question: An import of an Agilent array dataset from caArray into calIntegrator 2 failed. Why?

Topic: calIntegrator Usage

Release: calIntegrator 1.2 and above

Date entered: 08/31/2010

Answer

For Agilent array data, the probe names have to start with "A_". This lets the code ignore all of the other irrelevant data in a real Agilent design file.

Have a comment?

Please leave your comment in the [calIntegrator End User Forum](#).