

caNanoLab 2.3.11 Release Note

Welcome to the caNanoLab 2.3.11 Release Note, May 21, 2021. Topics in this Release Note include the following.

- [Improvements](#)
- [Fixed Bugs](#)
- [Known Issues and Workarounds](#)
- [Release History](#)



JIRA

If you have access to the JIRA Tracking System, click a JIRA link to review additional development details.

Improvements

The following table lists and describes the improvement(s) in this caNanoLab release.

Description of the improvement(s)	JIRA issue link
Upgraded AngularJS version to 1.8.	CANANOLAB-675

Fixed Bugs

The following table lists and describes the bug(s) resolved in this caNanoLab release.

Description of the bug fix(es)	JIRA issue link
Link changes on the caNanoLab Home Page.	CANANOLAB-88
Unable to delete composition files.	CANANOLAB-559
Advanced search not returning results for Sample Name Criteria: angular.js error.	CANANOLAB-677
Unable to search for samples.	CANANOLAB-687
Unable to create new samples.	CANANOLAB-688
Protocol: Unable to create new protocol when uploading\attaching file.	CANANOLAB-691
Publication: Unable to create new publication nor update publication when attaching an uploaded file.	CANANOLAB-692
Sample Composition: When creating composition file, on the Summary Page, the Title and Download Link for uploaded file is missing.	CANANOLAB-693
Sample: General Info: Unable to export JSON or XML file.	CANANOLAB-694
Sample: Characterization: Unable to add findings in any characterization element because the column name was not selectable.	CANANOLAB-695
Sample: Composition File: Unable to delete composition file.	CANANOLAB-697
Composition: Functionalizing Entity: Unable to edit or delete inherent function.	CANANOLAB-699
Composition: Nanomaterial Entity: Unable to edit or delete a composing element.	CANANOLAB-700

Characterization: Unable to create a characterization element with a technique and instrument added.	CANANOLAB-702
Unable to edit or delete publication authors.	CANANOLAB-704
Publications: Unable to edit or delete users or groups under Access to the Publication.	CANANOLAB-705
Samples: Unable to edit or delete users or groups.	CANANOLAB-706
Samples: Unable to delete existing samples.	CANANOLAB-707
Admin: Unable to add a role to a new user nor edit a user's existing role.	CANANOLAB-709
NCI-INC0572366-caNano account issues. When creating a new account, unable to reset password and change user role.	CANANOLAB-710

Known Issues and Workarounds

The following table lists and describes how to work around the following remaining issues in caNanoLab.

Known Issue	Workaround	JIRA issue link
The Results tab is populated when a curator/administrator submits a batch job to create a data availability matrix, but the tab disappears if the curator/administrator logs out and logs in again, even if the job is still running.	None	CANANOLAB-81
When submitting a new publication, clicking Reset does not reset the Submit Publication Form.	None	CANANOLAB-197
You cannot sort Protocol search results by the Protocol Type and Protocol Name columns.	None	CANANOLAB-367
You have to click the DOI external URL twice (or open the link in a new window by right-clicking) to view the external PubMed content.	None	CANANOLAB-368
A Publication Search returns more records than is valid based on the search criteria.	None	CANANOLAB-391
Export XML: Adding Other Ex Vv with special characters in the Characterization name prevents the Export xml from generating.	None	CANANOLAB-615
Export JSON: Special characters such as apostrophe and ampersand are showing up in the JSON export in their Unicode form.	None	CANANOLAB-616
Apostrophe does not appear correctly.	None	CANANOLAB-631
Composition: Unable to Edit nor Delete File from Nanomaterial Entity.	None	CANANOLAB-701
Characterization: Unable to Edit nor Delete a File from a Characterization Element's Finding, Can't Delete Findings with an attached File having Keywords. Root cause of issue is due to Adding Keywords when Uploading a File under Findings.	None	CANANOLAB-703
Sample Characterization: When creating new Model, Other Ex Vv, or Other, the Characterization Name is cleared after hitting the Save Button on Technique and Instruments and Findings.	None	CANANOLAB-708

Allow password reset without knowing old password.	Contact caNanoLab-Support@ISB-CGC.org .	CANAN OLAB-712
Characterization: General Page: Clean Up format for added Files In Findings: Missing Description Header.	None	CANAN OLAB-713
Composition: All Tab: Clean Up format for added Files In Findings: Bad spacing and Missing Headers.	None	CANAN OLAB-714
My Favorites: Deleted Items should be removed from My Favorites page.	None	CANAN OLAB-715

Release History

Release	Date
2.3.11	21 May 2021
2.3.10	30 October 2020
2.3.8	19 September 2019
2.3.7	11 June 2019
2.3.6	7 January 2019
2.3.5.1	15 November 2018
2.3.5	2 October 2018
2.3.4	13 July 2018
2.3.3	25 January 2018
2.3.2	October 2017
2.3	April 2017
2.2	January 2017
2.1.1	February 2016
2.1	November 2015
2.0	December 2014
1.5.5	May 2013
1.5.4	August 2012
1.5.3	May 2012
1.5.2	October 2010