

Defining Functionalizing Entity Properties

If you selected one of the Functionalizing Entity Types listed in the following table, the **Properties** section opens, and you must fill in additional information.

Functionalizing Entity Type	Properties to Complete
Antibody	<ol style="list-style-type: none">1. Select the Type of antibody functioning for this sample.<ul style="list-style-type: none">• Fab• ScFv• whole2. Select the Isotype for this antibody.<ul style="list-style-type: none">• IgA• IgD• IgE• IgG• IgM3. Select the animal Species from which the antibody was derived.
Biopolymer	<ol style="list-style-type: none">1. Select the Biopolymer type of antibody functionalizing entity for this sample. 2. Using the appropriate format, enter the complete Sequence of the biopolymer.
Small Molecule	Add an Alternate Name for the Small Molecule.