

Menu Options

The following table lists and describes each caNanoLab menu.

Menu Options	Description
Home	Returns to the home page.
Workflow	Presents a workflow graphic with active links to access caNanoLab functions. For details, refer to Using the Workflow to Get Started
Protocols	Search for existing protocols, create a new protocol, and submit a protocol for review for public access. For details, refer to Managing Protocols in caNanoLab .
Samples	<p>Search for existing samples, create a new sample, and submit a sample for review for public access, as described in Managing Samples in caNanoLab. When you add a sample, the following options appear in the Navigation Tree:</p> <div data-bbox="253 558 479 785"> <div>NAVIGATION TREE</div> <div>GENERAL INFO</div> <div>COMPOSITION</div> <div>CHARACTERIZATION</div> <div>PUBLICATION</div> </div> <ul style="list-style-type: none"> • Composition to add annotations and publications to the sample. For details, refer to Managing Sample Composition Annotations. • Characterization to review or ascribe characterizations to samples. For details, refer to Managing Sample Characterizations. • Publication to add publication to a sample. For details, refer to Adding Publications to a Sample.
Publications	Search for sample publications, create a new publication, or submit a publication for review for public access. For details, refer to Managing Publications in caNanoLab .
Groups	<p>Create and manage collaboration groups which allow researchers to share data among members of the group.</p> <ul style="list-style-type: none"> • <i>Curators</i> can manage all collaboration groups and view a list of samples assigned to each group. • If you do not have the curator role, you can manage only your own collaboration groups. <p>For details, refer to Managing Collaboration Groups.</p>
Curation	<i>Curator only:</i> Review items pending public release and generate sample data metrics in a batch. For details, refer to Managing Data Curation .
Results	<i>Curator only:</i> Appears after you perform a batch process. For details, refer to Reviewing Batch Results .
My Workspace	Review the samples, protocols, and publications that you saved or are associated with your log in. For details, refer to Using My Workspace .
My Favorites	Review the samples, protocols, and publications that you bookmarked for easy access. For details, refer to Using My Favorites .
Admin	<i>Admin only:</i> Add and edit users. For details, refer to Managing caNanoLab User Accounts .
Login /Logout	Login or out of caNanoLab.