Searching for a Template from AIM Template Service 1.0

You can search the AIM Template Service for a template that you can customize in ATB. Once you find the template you want to use in the search results, you can import it.

1. Click the Template tab.

? Unknown Attachment

and then select Search Templates.

2. Below the left panel, click
The Template Search page appears.

? Unknown Attachment

3. Narrow your search by entering information into the fields. If you do not enter any information into the fields, all template groups available in the AIM Template Service appear, and it may take longer to see the results of your search. The following fields are available.

Field	Description
Template Name	Template name
Template Code Desc	Template code description
Template Description	Template description
Template Author	Person who created the template in the AIM-enabled application
Disease	Disease that is the topic of the templates in the template group
Body Part	Body part that is the topic of the templates in the template group
Modality	Any of the various types of equipment or probes used to acquire images of the body, such as radiography, ultrasound and magnetic resonance imaging

4. Click Search.

All templates matching the search criteria you entered appear.

? Unknown Attachment