Adding a Template Group 2.0

1. Click the Template Groups tab.

	elow the left panel, clic	
Th	e Group Members panel at the bottom of the application expands to include a new template group panel.	
and the second	Name* :	
and the second	Version* :	
	Description* :	
	Authors* :	Default Author
and the second	UID* :	2.25.131995507397654214847534467268280544622
and the second	Date :	2013-01-14
Save		

3. Enter the name of the new group, its version number and description, and your name in the boxes provided. You can edit this information later if needed. 4. Click **Save**.

The new template group appears in the Groups panel on the left side of the application.