



Creating a Characteristic Quantification 2.0

1. [Add a template component](#) of the Imaging Observation or Anatomic Entity type.
2. [Add a characteristic to a template component](#) of the Imaging Observation or Anatomic Entity type.
3. [Add an answer choice](#) to the characteristic.
4. Right-click the answer choice and select **Create Characteristic Quantification** ( Unknown Attachment).
5. In the panel on the right, enter a name for the characteristic and select [Annotator Confidence](#), if that attribute applies to the template component.
6. Click **Save**.
7. Create a characteristic quantification type. The options include [Non Quantifiable](#), [Scale](#), [Numerical](#), [Quantile](#), and [Interval](#). Except for the Non Quantifiable type, which can always be added, you cannot create more than one type for a single answer choice.

Note

Delete any quantification from an answer choice by right-clicking it and selecting

 Unknown Attachment